PCT

WORLD INTELLECTUAL PROPERTY ORGANIZATION International Bureau



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁶:

A61L 2/00

A2

(11) International Publication Number: WO 98/30249

(43) International Publication Date: 16 July 1998 (16.07.98)

(21) International Application Number: PCT/US98/00635

(22) International Filing Date: 13 January 1998 (13.01.98)

(30) Priority Data: 08/782.647 14 January 1997

14 January 1997 (14.01.97) US

(71) Applicant: STERIS CORPORATION [US/US]; 5960 Heisley Road, Mentor, OH 44060 (US).


- (72) Inventors: HAIRE, Donna, M.; 5273 E. 117th Street, Garfield Heights, OH 44125 (US). WOLF-BLOOM, David, F.; Apartment 203, 2540 North Moreland Boulevard, Shaker Heights, OH 44120 (US). MALCHESKY, Paul, S.; 239 Barrington Ridge, Painesville Twp., OH 44077 (US).
- (74) Agent: KOCOVSKY, Thomas, E., Jr.; Fay, Sharpe, Beall, Fagan, Minnich & McKee, Suite 700, 1100 Superior Avenue, Cleveland, OH 44114–2518 (US).

(81) Designated States: AU, CA, JP, KR, SG, European patent (AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE).

Published

Without international search report and to be republished upon receipt of that report.

(54) Title: METHOD AND APPARATUS FOR CLEANING, DECONTAMINATING, AND STERILIZING CATHETERS



(57) Abstract

A method for cleaning, decontaminating, and sterilizing catheters (10) using a combination of liquid and gaseous/plasma sterilization techniques to ensure the complete and efficient sterilization of a catheter (10). Angiographic dye and saline are removed from the interior (36) of the balloon (14) and its lumen (16). The outer surfaces of the catheter (10) and a guide wire lumen (18) of the catheter (10) are cleaned, decontaminated, and sterilized (42) with a liquid sterilant. The liquid sterilant fills a balloon (14) and a balloon lumen (16) of the catheter (10). The liquid sterilant is retained in the balloon (14) and the lumen (16) for a selected amount of time. Thereafter, the liquid sterilant is drained from the balloon (14) and the balloon lumen (16). The filling, retaining, and draining steps are repeated until an interior (36) of the balloon (14) and the balloon lumen (16) are sterilized. Residual liquid sterilant is rinsed from the interior (36) of the balloon (14) and the balloon lumen (16). The catheter (10) is dried and then a plasma or gaseous sterilant is used to sterilize at least the outer surfaces and the guide wire lumen (18) of the catheter (10).

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

| AL | Albania | ES | Spain | LS | Lesotho | SI | Slovenia |
|---------------|--------------------------|----|---------------------|------------------------|-----------------------|----|--------------------------|
| AM | Armenia | FI | Finland | LT | Lithuania | SK | Slovakia |
| AT | Austria | FR | France | LU | Luxembourg | SN | Senegal |
| AU | Australia | GA | Gabon | LV | Latvia | SZ | Swaziland |
| AZ | Azerbaijan | GB | United Kingdom | MC | Monaco | TD | Chad |
| BA | Bosnia and Herzegovina | GE | Georgia | MD | Republic of Moldova | TG | Togo |
| BB | Barbados | GH | Ghana | MG | Madagascar | TJ | Tajikistan |
| BE | Belgium | GN | Guinea | MK | The former Yugoslav | TM | Turkmenistan |
| \mathbf{BF} | Burkina Faso | GR | Greece | | Republic of Macedonia | TR | Turkey |
| BG | Bulgaria | HU | Hungary | ML | Mali | TT | Trinidad and Tobago |
| ВJ | Benin | IE | Ireland | MN | Mongolia | UA | Ukraine |
| BR | Brazil | IL | Israel | MR | Mauritania | UG | Uganda |
| BY | Belarus | IS | Iceland | MW | Malawi | US | United States of America |
| CA | Canada | IT | Italy | MX | Mexico | UZ | Uzbekistan |
| CF | Central African Republic | JP | Japan | NE | Niger | VN | Viet Nam |
| \mathbf{CG} | Congo | KE | Kenya | NL | Netherlands | YU | Yugoslavia |
| CH | Switzerland | KG | Kyrgyzstan | NO | Norway | zw | Zimbabwe |
| CI | Côte d'Ivoire | KP | Democratic People's | NZ | New Zealand | | |
| CM | Cameroon | | Republic of Korea | PL | Poland | | |
| CN | China | KR | Republic of Korea | PT | Portugal | | |
| CU | Cuba | KZ | Kazakstan | RO | Romania | | |
| CZ | Czech Republic | LC | Saint Lucia | RU | Russian Federation | | |
| DE | Germany | LI | Liechtenstein | SD | Sudan | | |
| DK | Denmark | LK | Sri Lanka | SE | Sweden | | |
| EE | Estonia | LR | Liberia | $\mathbf{s}\mathbf{G}$ | Singapore | | |
| | | | | | | | |
| | | | | | | | |

- 1 -

Method and Apparatus for Cleaning, Decontaminating, and Sterilizing Catheters

Background of the Invention

The present invention relates to the art of cleaning, decontaminating, and sterilizing medical instruments and devices. It finds particular application in conjunction with the cleaning, decontamination, and sterilization of previously used catheters, such as percutaneous transluminal coronary angioplasty catheters (PTCA catheters), for subsequent reuse and will be described with particular reference thereto.

10 As the costs associated with the delivery of effective healthcare continue to rise, there is a continuing effort on the part of doctors, healthcare facilities, and insurers to find ways of reducing these costs. The increasing popularity of single-use 15 disposable medical devices has contributed to the increasing costs of surgery. The PTCA catheter is a widely used single-use disposable medical device used by surgeons to perform percutaneous transluminal angioplasty, a surgical procedure to clear blockages of 20 plaque and the like from the arteries of a patient.

A percutaneous transluminal coronary angioplasty is carried out utilizing a PTCA catheter which typically takes the form of a balloon and a balloon lumen having a coaxial guide wire lumen passing therethrough. The proximal end of the PTCA catheter is bifurcated and includes a guide wire lumen port providing access to the guide wire lumen and a balloon lumen port providing access to the balloon lumen, and thus, the balloon. The distal end of the PTCA catheter includes the closed balloon which is sealed around the open end of

- 2 -

the guide wire lumen, i.e., the guide wire lumen passes through the balloon and is open at its distal end. The surgeon uses a guide wire to guide the PTCA catheter into the vascular system of the patient to an area where a blockage exists. The catheter is often inserted into the femoral artery of the patient, and therefore, the catheter must be of a significant length (often 150 centimeters) to reach areas near the patient's heart.

After the catheter is in position, the surgeon inflates the balloon with a fluid, such as saline, by injecting the saline into the balloon lumen port of the catheter. Often, the saline includes a radiographic dye which enables the position and the inflated profile of the balloon to be monitored radiographically. The balloon expands radially and thus forces the arterial plaque radially outward into the wall of the artery. The surgeon then deflates the balloon, withdraws the catheter from the patient's vascular system, and disposes of the catheter.

Usually, more than one catheter is used on a particular patient, at a cost of several hundred dollars or more per catheter. Thus, a significant cost savings could be realized from the reuse of PTCA catheters. One difficulty associated with cleaning and sterilizing catheters such as PTCA catheters is attributable to a combination of the length of the catheter, the narrowness of the catheter and the guide wire and balloon lumens, and the closed balloon at the distal end of the catheter which does not allow liquid sterilant to be flushed through the balloon lumen and the balloon.

The cleaning of these devices by hand is time consuming and unreliable. Prior attempts to automate the process have resulted in unduly complex systems that have

- 3 -

not been proven to be effective. Prior systems rely solely on liquid sterilization and are not compatible with gas sterilization techniques such as the use of ethylene oxide gas (EO) and these prior systems 5 specifically teach away from the use of EO. Even with respect to the liquid sterilization utilized in the prior systems, the methods and apparatus disclosed do not control the "contact time" of the cleaning agent and the sterilant within the balloon and balloon lumen. They 10 require the balloon to be pressurized with liquid sterilant and a cleaning agent an excessive number of times, reducing the remaining useful life of the balloon. Furthermore, merely filling the balloon lumen and balloon with a cleaning agent and sterilant and immediately 15 flushing the same therefrom a selected amount of times, without further controls, has been found to be an ineffective cleaning and sterilization technique in certain circumstances.

One prior PTCA catheter sterilization method requires the liquid sterilized catheter to be rinsed of liquid sterilant immediately prior to use. This prior method requires a nurse or a surgical assistant to rinse the balloon and balloon lumen by injecting several syringes or more of sterile saline into the balloon. The saline that drains from the balloon and balloon lumen must be tested for the presence of sterilant using chemical test strips. The rinsing must be continued in this manner until the test strips indicate that the concentration of sterilant in the evacuated rinse fluid is below safe thresholds. It can be seen that the foregoing PTCA catheter sterilization method is time consuming and inconvenient. Also, the rinsing of the catheter just prior to its use results in a wet catheter

- 4 -

being presented to the surgeon, which some surgeons may find objectionable.

Another disadvantage of known catheter sterilization systems is their requirement that a single 5 catheter be processed or sterilized and stored in a specially designed cassette while awaiting use. These cassettes must be specially pressurized and capped in order for them to maintain the sterility of the catheter. Also, each cassette holds only one catheter, and thus, 10 the efficiency of this prior system is necessarily limited.

The present invention provides a new and improved catheter cleaning, decontamination, and sterilization method and apparatus which overcomes the above-referenced problems and others.

Summary of the Invention

In accordance with one aspect of the present invention, a method of cleaning, decontaminating, and sterilizing an elongated catheter having outer surfaces, a balloon, and a dead-end balloon lumen open at a first end and closed by the balloon at a second end is provided. Using a liquid decontaminant, the outer surfaces of the catheter are cleaned and decontaminated.

- 25 The balloon lumen and the interior of the balloon are cleaned, decontaminated, and sterilized with a liquid sterilant. The catheter is dried and thereafter, at least the outer surfaces of the catheter are sterilized with a gaseous or plasma sterilant.
- In accordance with a more limited aspect of the present invention, the step of sterilizing the balloon lumen and the interior of the balloon of the catheter includes rinsing the dye and saline from within the

- 5 -

balloon and the balloon lumen. Liquid sterilant is communicated into the interior of the balloon and into the balloon lumen and retained in the balloon and balloon lumen for a selected amount of time. The liquid sterilant is then drained from the balloon and balloon lumen of the catheter. The filling, retention, and draining steps are repeated until the interior of the balloon and the balloon lumen are sterilized.

In accordance with yet another aspect of the 10 present invention, an apparatus for cleaning, decontaminating, and sterilizing at least the balloon lumen and the interior of the balloon of a plurality of catheters using the method of the present invention is provided. The apparatus includes a supply of liquid 15 sterilant and a liquid bath for receiving a plurality of catheters therein to maintain the catheters at a selected temperature. A manifold is in selective fluid communication with the liquid sterilant supply and includes a plurality of catheter connection fittings for 20 connecting the balloon lumen of each of the plurality of catheters in fluid communication with the manifold. pump pumps liquid sterilant under pressure from the supply of liquid sterilant into the balloon lumen and balloon of each of the plurality of catheters connected 25 to the manifold and also evacuates the liquid sterilant from the balloon and the balloon lumen of each of the plurality of catheters. At least one valve selectively blocks the flow of liquid sterilant to and from the manifold. An electronic controller is connected to the 30 pump and the at least one valve for controlling the operation of the pump and the at least one valve. controller is programmed to control the pump and the at least one valve to repeat the operations of operating the

- 6 -

pump to fill the balloon and balloon lumen of each of the plurality of catheters connected to the manifold with liquid sterilant; closing the at least one valve to retain the liquid sterilant within the balloon lumen and balloon of each of the plurality of catheters connected to the manifold for a selected amount of time; and, opening the at least one valve and operating the pump to evacuate essentially all of the liquid sterilant from within the balloon and balloon lumen of each of the plurality of catheters connected to the manifold.

One advantage of the present invention is the more effective cleaning and sterilization of catheters.

Another advantage of the present invention is the efficiency and ease with which catheters are sterilized.

Another advantage of the present invention is the combination of liquid and gas or plasma sterilization techniques to effectively sterilize catheters.

Another advantage of the present invention is the 20 elimination of the need to sterilize and store the catheters in special pressurized containers.

Another advantage of the present invention is that the sterilized catheters are rinsed, dried, and ready for immediate use without further rinsing or 25 preparation.

Another advantage of the present invention is the ability to simultaneously sterilize the interior of the balloon and the balloon lumen of each of a plurality of catheters.

30 Still further advantages of the present invention will become apparent to those of ordinary skill in the art upon reading and understanding the following detailed description of the preferred embodiments.

- 7 -

Brief Description of the Drawings

The invention may take form in various components and arrangements of components, and in various 5 steps and arrangements of steps. The drawings are only for purposes of illustrating a preferred embodiment and are not to be construed as limiting the invention.

FIGURE 1A shows a PTCA catheter suitable for being sterilized in accordance with the method and 10 apparatus of the present invention;

FIGURE 1B is a partial enlarged diagrammatic illustration of distal portion of the PTCA catheter of FIG. 1, including the balloon, balloon lumen, and guide wire lumen thereof;

FIGURE 2 is a flow chart illustrating catheter cleaning and sterilization in accordance with the present invention;

 $\,$ FIGURE 3 is a flow chart illustrating removing the catheter from the patient and preserving the catheter

20

for sterilization in accordance with the present invention:

FIGURE 4 is a flow chart illustrating performing the first cleaning and decontamination of the 25 catheter:

FIGURE 5 is a flow chart illustrating cleaning, decontaminating, and sterilizing the balloon lumen and interior of the balloon;

FIGURE 6 is a flow chart illustrating the 30 gaseous or plasma sterilization of the catheter;

FIGURE 7 is a schematic view of a balloon and balloon lumen pump system as it is configured for liquid sterilization of the balloon and balloon lumen; and

- 8 -

FIGURE 8 is a schematic view of a balloon and balloon lumen pump system as it is configured for rinsing the balloon and balloon lumen.

5 <u>Detailed Description of the Preferred Embodiments</u>

With reference to FIGURES 1A and 1B, a percutaneous transluminal coronary angioplasty (PTCA) catheter 10 includes an elongated and flexible tubular member 12 with a balloon 14 located at the distal end.

- 10 As is seen most clearly in FIGURE 1A, the tubular member 12 provides a channel or lumen 16 to the balloon 14. A through-lumen, referred to as a guide wire lumen 18, is coaxially positioned within balloon lumen 16 and extends through balloon 14 to an open distal end 20. The balloon
- 15 14 is closed at its distal end 22 around the guide wire lumen 18.

The tubular member 12 of the catheter 10 has a proximal end 24 including a bifurcated hub 26. bifurcated hub 26 includes a guide wire lumen access port 20 28 providing access to the guide wire lumen 18 and a balloon lumen access port 30 providing access to the balloon lumen 16, and thus to the balloon 14. Each port 30 typically includes commonly known luer lock connectors 32 or the equivalent. As shown herein, the 25 bifurcated hub 26 is connected to the tubular member 12 through a strain relief sheath 34 for a more durable interconnection. The illustrated PTCA catheter 10 is an example of a catheter that can be cleaned, decontaminated, and sterilized utilizing a method and 30 apparatus in accordance with the present invention.

Those skilled in the art will recognize that the PTCA catheter 10 may be modified but will remain suitable for

- 9 -

sterilization in accordance with the present invention, and that a sterilization method and apparatus in accordance with the present invention may be utilized to clean and sterilize other catheters having a balloon and balloon lumen as described above.

During angioplasty procedures, the tubular member 12 of the catheter 10 is inserted into the vascular system of the patient, and thus, the outer surfaces of the catheter 10, including the outer surfaces 10 of the hub 26, the strain relief sheath 34, the tubular member 12, and the balloon 14 become contaminated with the blood and other bodily fluids of the angioplasty patient. Furthermore, the guide wire lumen 18, which is open at the distal end 20, also becomes contaminated with 15 the patient's blood. In most angioplasty procedures, the balloon lumen 16 and the interior 36 of the balloon 14 will not become contaminated with the patient's bodily fluids. In subsequent angioplasty procedures, any contamination within the balloon 14 would not be readily 20 transmissible into the blood of the later patient in light of the closed nature of the balloon. However, the balloon 14 may burst upon inflation in a later procedure, and for this and other reasons, the safe sterilization of PTCA catheters 10 includes a complete sterilization of 25 the balloon lumen 16 and the interior 36 of the balloon 14.

Furthermore, during angioplasty procedures, an angiographic dye and saline mixture is communicated into the interior 36 of the balloon 14. This dye and saline 30 mixture are removed from the interior 36 of the balloon 14 within a relatively short time (less than 24 hours) lest it crystallizes, adversely affecting the utility of

- 10 -

the catheter 10 for subsequent angioplasty procedures. The method and apparatus of the present invention result in the sterilization of an entire PTCA catheter 10, including all of the outer surfaces, as well as the guide 5 wire lumen 18, the balloon lumen 16, and the interior 36 of the balloon 14.

With reference to FIGURE 2, the preferred overall method of sterilizing a PTCA catheter 10 is shown. Each step or means shown in FIGURE 2 is described 10 in further detail with reference to FIGURES 3 - 6. step 40 the catheter 10 is recovered from the patient and preserved for sterilization. Once the catheter 10 has been recovered and preserved for later sterilization, a step or means 42 cleans and decontaminates the outer 15 surfaces and the guide wire lumen 18 of the catheter 10. Following the cleaning and decontamination of the outer surfaces and the guide wire lumen 18 of the catheter 10, a step or means 44 cleans, decontaminates, and sterilizes the balloon lumen 16 and the interior 36 of the balloon 20 14 of the catheter 10. Preferably, a step or means 46 provides a second cleaning and decontamination of the outer surfaces and the guide wire lumen 18 of catheter 10. A step or means 48 dries the catheter 10, and a step or means 50 sterilizes the catheter with a 25 gaseous sterilant, resulting in a completely sterilized catheter that is free of harmful amounts of residual sterilant, dry, and ready for reuse.

With reference to FIGURE 3, the step 40 of recovering the catheter 10 from the patient is shown in 30 greater detail and includes steps or means 40a - 40d. In a step 40a, the catheter 10 is removed from the patient. This step is generally performed manually. Immediately

- 11 -

following the removal of the catheter 10 from the patient, a step or means 40b fills the balloon 14 of the catheter 10 with a sterile liquid, such as sterile saline, sterile water, or the like, and caps the balloon 1 lumen port 30 of the catheter 10 to retain the sterile liquid in the balloon 14 and balloon lumen 16. The sterile liquid in the balloon 14 and balloon lumen 16 and the capping of the lumen port 30 prevents the angiographic dye and saline within the balloon 14 from 10 crystallizing. Furthermore, capping the balloon lumen port 30 not only prevents the escape of the sterile fluid from the balloon 14, but also prevents the contamination of the interior 36 of the balloon 14 with blood and other contaminates that are in the guide wire lumen 18 and on 15 the outer surfaces of the catheter 10.

A step or means 40c flushes the guide wire lumen 18 with a liquid and a step or means 40d submerges or otherwise surrounds the catheter 10 with a liquid. Although ordinary saline or the like may be utilized as 20 the liquid in the step or means 40c and the step or means 40d, the liquid utilized is preferably an anticoagulant liquid solution. One suitable solution used in the steps or means 40c and 40d is a solution of 10,000 units of Heparin per 1000 milliliters (ml.) of saline. With such 25 a solution, the guide wire lumen 18 of the catheter 10 is flushed with approximately 10 ml. of the solution and the catheter is then submerged in the solution for about approximately 15 minutes. The anticoagulant solution rinses blood and the like from the outer surfaces of the 30 catheter 10 and from within the guide wire lumen 18. also prevents the clotting and hardening of the blood on the catheter 10 and in the guide wire lumen 18.

- 12 -

skilled in the art will recognize that the capped balloon lumen port 30 prevents the contamination of the interior 36 of the balloon 14 with blood when the catheter 10 is treated with the anticoagulant solution. Once the catheter 10 has been recovered from the patient and treated as described, an indefinite amount of time may pass. The catheter 10 can remain submerged in a bath of liquid as described. In practice, the catheter sterilization process continues relatively quickly following the recovery of the catheter 10 from the patient.

With reference to FIGURE 4, the sterilization process continues with the step or means 42 cleaning and decontaminating the outer surfaces and the guide wire 15 lumen 18 of the catheter 10. The balloon lumen port 30 remains capped at this point in the sterilization process. A step or means 42a removes the catheter 10 from the anticoagulant liquid and a step or means 42b surrounds the catheter 10 with liquid sterilant, decontaminant, or disinfectant and flushes the guide wire lumen 18 with liquid sterilant, decontaminant, or disinfectant for a selected amount of time.

Thereafter, a step or means 42c rinses the outer surfaces and the guide wire lumen 18 of the catheter 10 with a sterile rinse liquid such as sterile water. While the step or means 42 decontaminates the outer surfaces of catheter 10 and also facilitates the decontamination of the guide wire lumen 18, additional steps are provided as described below to ensure the decontamination of these areas of catheter 10. A preferred apparatus for rinsing and decontaminating the interior surfaces of the catheter 10 is illustrated in

U.S. Patent Nos. 4,731,222, 4,892,706, 5,116,575, and 5,217,698 of an assignee herein.

The process either continues immediately with the step or means 44 to sterilize the balloon lumen 16 and the interior 36 of the balloon 14 or a step or means 42d places the catheter 10 in a suitable temporary storage container, such as a sealable biohazard bag, for later processing. The step or means 42d is utilized when further processing is conducted at a different site.

10 Also, a hospital may wish to accumulate a large number of catheters for subsequent processing in a batch as is described below.

As is shown in detail in FIGURE 5, reprocessing continues with the step or means 44 to 15 clean, decontaminate, and sterilize the balloon lumen 16 and the interior 36 of the balloon 14. A step or means 44a uncaps the balloon lumen port 30 and withdraws the sterile fluid and any air from the balloon 14 and balloon lumen 16. Because the angiographic dye and saline 20 communicated into the interior 36 of the balloon 14 and the balloon lumen 16 during the angioplasty procedure can interfere with the biocidal action of the liquid sterilant and the further use of the catheter 10, the sterilization of the balloon 14 and the balloon lumen 16 25 involves the removal of the dye and saline. A step or means 44b fills the balloon 14 and the balloon lumen 16 with a sterile rinse liquid, such as sterile water. sterile water may optionally be retained within the balloon 14 and the lumen 16 for a selected "contact time" 30 to facilitate the rinsing process. A step or means 44c evacuates the sterile rinse fluid from the balloon 14 and the balloon lumen 16 and a step or means 44d repeats the

- 14 -

filling and evacuating steps 44b, 44c until the dye and saline are removed from the balloon 14 and the balloon lumen 16. It has been found that repeating the steps 44b, 44c at least fifteen and preferably twenty times, with an optional contact time of approximately fifteen - thirty seconds is effective for removing the dye and saline from the interior 36 of the balloon 14 and from the balloon lumen 16.

In order to sterilize the interior 36 of the 10 balloon 14 and the balloon lumen 16, a liquid sterilant is fed into the balloon 14 and the balloon lumen 16. While simply repeatedly filling the balloon 14 and the lumen 16 with liquid sterilant and immediately evacuating the same may sterilize the interior 36 of the balloon 14 15 and the lumen 16, it has been found that such a method does not sterilize the interior 36 of the balloon 14 and the balloon lumen 16 as effectively as repeatedly filling the balloon 14 and balloon lumen 16 with the liquid sterilant, retaining the liquid sterilant in the balloon 20 14 and balloon lumen 16 for a selected amount of time, and thereafter evacuating the liquid sterilant from the balloon 14 and the balloon lumen 16. This increase in "contact time" or the amount of time that the liquid sterilant is in contact with the interior 36 of the 25 balloon 14 and the balloon lumen 16 results in more effective sterilization. It has also been found that filling the balloon 14 once with liquid sterilant and holding it filled for a selected duration is less effective than filling and emptying the balloon several 30 times over the same selected duration.

Therefore, a step or means 44e fills the

- 15 -

balloon 14 and the balloon lumen 16 with a liquid sterilant and maintains the filled state of the balloon 14 for a selected amount of time. A step or means 44f evacuates the liquid sterilant from the balloon 14 and 5 lumen 16 and a step or means 44g repeats the step 44e and the step 44f until the interior 36 of the balloon 14 and the balloon lumen 16 are sterilized. For example, it has been found that 106 colony forming units of spores are killed or removed from the interior 36 of the balloon 14 10 and the balloon lumen 16 by repeating the steps 44e and 44f at least eight times with a contact time of at least approximately two minutes. It is preferable to repeat the steps 44e and 44f at least ten times using a contact time of at least approximately one and one half minutes 15 using a 0.2% peracetic acid sterilant. By way of example, the interior 36 of a balloon 14 and a balloon lumen **16** of a catheter **10** were inoculated with $10^6 - 10^7$ / 0.5 ml. B. subtilis spores using different carrier solutions of water, saline, and 10% serum / 90% saline. 20 A 0.2% peracetic acid sterilant utilized. was Sterilization was not achieved, with the water carrier, when steps 44e and 44f were repeated as follows: one time with a contact time of five minutes; two times with a contact time of seven minutes; three times with a contact 25 time of one minute; five times with a contact time of one minute; and, ten times with a contact time of 1 minute (3 out of four trials not sterile). Sterilization was achieved, with the water carrier, when steps 44e and 44f were repeated as follows: four times with a contact time 30 of two minutes; twenty-two times with a contact time of 0.25 minutes; and, twenty-four times with a contact time of 0.25 minutes. Using the serum and saline carrier

solution which presents the most difficult sterilization situation, when steps 44e and 44f were repeated ten times with a contact time of at least 1.5 - 2 minutes, sterilization resulted for seven of the nine trials.

5 With the saline carrier, six trials were conducted and steps 44e and 44f were repeated ten times with a contact time of two minutes. All six trials resulted in sterilization.

Following the sterilization of the interior 36 10 of the balloon 14 and the balloon lumen 16 using the liquid sterilant as described, the liquid sterilant is rinsed from the interior 36 of the balloon 14 and from the balloon lumen 16. Rinsing the interior 36 of the balloon 14 and the balloon lumen 16 is preferably 15 accomplished by repeatedly filling the balloon 14 and the balloon lumen 16 with a sterile rinse liquid, such as sterile water, and thereafter immediately evacuating the rinse liquid from the balloon 14 and the balloon lumen A step or means 44h fills the balloon 14 and the 20 balloon lumen 16 with a sterile rinse liquid. A step or means 44i evacuates the sterile rinse liquid from the balloon 14 and lumen 16 after the selected amount of time has passed. A step or means 44j causes the means 44h and 44i to repeat the filling and evacuating steps until the 25 amount of residual sterilant in the interior 36 of the balloon 14 and in the balloon lumen 16 is below the selected threshold.

For example, using a peracetic acid sterilant and a sterile rinse liquid of sterile deionized/filtered 30 water, it has been found that repeating steps 44h and 44i at least fifteen times, and preferably twenty times, results in the interior 36 of the balloon 14 and the

lumen 16 being essentially free of residual sterilant.

Once the rinsing process is complete, a step or means 44k evacuates the sterile rinse liquid from the balloon 14 and its lumen 16 and a step or means 44l caps the balloon 5 lumen port 30 to prevent recontamination.

In case of any contamination of the catheter 10 during transport or handling with respect to the above procedures, it is preferable to provide an optional step or means 46 for a second cleaning and decontamination of 10 the outer surfaces and the guide wire lumen 18 of the catheter 10. The step or means 46 performs substantially the same operations as the step or means 42 as described above in detail with reference to FIGURE 4. The step or means 46 ensures that the outer surfaces of the catheter 15 10 are decontaminated and facilitates the sterilization of the guide wire lumen 18. Following the second cleaning and decontamination of the outer surfaces and the guide wire lumen 18 of the catheter 10 with a liquid sterilant, the catheter 10 is completely dried by a step 20 or means 48 in preparation for a step or means 50 to gas sterilize the catheter 10. As is shown in FIGURE 6, a step or means 48a uncaps the balloon lumen port 30 of the catheter 10. A step or means 48b dries the outer surfaces, the guide wire lumen 18, the balloon lumen 16, 25 and the interior 36 of the balloon 14. The step or means 48b is preferably carried out in a sterile vacuum oven system with -25 to -30 inches of mercury (Hg) at a temperature of 30 degrees C to 40 degrees C for 1.5 hours. Alternatively, the oven can be a non-sterile 30 vacuum oven system.

With continuing reference to FIGURE 6, once the

catheter 10 is dry, the step or means 50 gas sterilizes the catheter 10, and in particular the outer surfaces and the guide wire lumen 18 thereof, preferably using a gaseous or plasma sterilant such as ethylene oxide (EO) 5 gas, vapor phase hydrogen peroxide, or the like. A step or means 50a deflates the balloon 14 and places a protective sheath over the balloon 14. A step or means 50b places the catheter in gas sterilization packaging, such as a standard EO sterilization bag. A step or means 10 50c exposes the packaged catheter 10 to a plasma or gaseous sterilant such as EO for a selected amount of time. For example, the step or means 50c exposes the catheter 10 to 100% EO for approximately 4 to 5 hours at 90 degrees F. A 4.75 hour EO exposure time has been 15 found to be effective. Following the exposure of the catheter 10 to EO, a step or means 50d aerates the catheter 10 and the packaging to disperse residual gaseous sterilant. It has been found effective for the step or means 50d to aerate the catheter 10 and the gas 20 sterilization packaging for 12 hours at approximately 90 degrees F. The catheter 10 can then be safely reused for another angioplasty procedure with a different patient.

With reference to FIGURE 7, a supply of liquid sterilant 62 is heated in a water bath 64. The heated 25 water bath 64 also receives a plurality of the catheters 10 during the sterilization process. The liquid sterilant, such as the 0.2% peracetic acid, is utilized at a temperature of 40 - 60 degrees C, preferably 50 - 56 degrees C. The heated water bath 64 heats the supply of liquid sterilant 62 and the catheters 10 placed in the bath 64 to the optimum operative temperature of the selected sterilant. A manifold 66 includes a plurality

of catheter connection fittings **68** for fluidically connecting the balloon lumen port **30** of each of the plurality of catheters **10** to the manifold **66**.

The manifold 66 is in fluid communication with 5 the supply of liquid sterilant 62 through a pump 70 and is also in fluid communication with a liquid sterilant disposal reservoir 72 through the pump 70. The liquid sterilant is pumped by the pump 70 from the supply 62 into the manifold 66 as indicated by the arrows \mathbf{A} and is 10 pumped from the manifold 66 to the disposal reservoir 72 as indicated by the arrows B. An inlet valve 74 and an outlet valve 76 selectively block the flow of fluid to and from manifold 66. A pressure transducer 78 monitors the pressure of the fluid in the manifold 66. A display 15 80 displays process parameters and information regarding the balloon and balloon lumen sterilization process to the operator of the apparatus 60. The parameters preferably include the date, the time, the sterilant and water bath temperatures, the number of balloon inflation 20 cycles, the duration of each balloon inflation cycle, and the like. Also, the apparatus is preferably enclosed within a fume hood or another suitable enclosure 86 to minimize the escape of sterilant vapors into the atmosphere, and to minimize environmental contamination 25 of the water bath 64, the supply of sterilant 62, and the supply of sterile rinse fluid 82 (FIGURE 8).

of the apparatus of FIGURE 7 during the rinsing steps 44b - 44d and 44h - 44j in which the angiographic dye and saline mixture, and the liquid sterilant, are respectively rinsed from the interior 36 of the balloons 14 and the balloon lumens 16 of each catheter 10

connected to the fittings 68 of manifold 66. A supply of sterile rinse liquid 82 is in fluid communication with the manifold 66 through the pump 70 and the manifold 66 is also in fluid communication with the disposal reservoir 72. The sterile rinse liquid is pumped from supply 82 to manifold 66 through conduits as shown by the arrows A, and the sterile rinse liquid is pumped from the manifold 66 to the disposal reservoir 72 through conduits as indicated by the arrows B. For rinsing operations, 10 each catheter 10 connected to the manifold 66 is preferably submerged in an ultrasonic water bath 84.

The water baths 64,84, the pump 70, the valves 74,76, the pressure transducer 78, and the visual display 80 are all electrically connected to an electronic 15 controller 88 that controls the operation of each of the noted components. The electronic controller 88 controls the application of heat and ultrasonic waves to the baths The electronic controller 88 also selectively operates the pump 70 as required to pump sterilant and 20 sterile rinse liquid, and monitors the fluid pressure in the manifold 66 through electrical signals received from the transducer 78. The controller 88 selectively opens and closes the valves 74,76 to selectively block the flow of liquid sterilant and rinse liquid to and from the 25 manifold 66 and the disposal reservoir 72, and displays process parameters and information to an operator on the display 80. An input system 90 is connected to the electronic controller 88 for operator input of data and for operator control. A printer 92 prints out the date 30 and time, the operator, and cycle information to provide written documentation for each catheter 10. The cycle

- 21 -

information includes such information as the number of balloon inflations with the liquid sterilant, the duration of sterilant contact, the number of rinse cycles, and the like.

5 The electronic controller 88 is programmed to operate the apparatus in accordance with the step 44 as described above with reference to FIGURES 2 and 5. Through the operation of the pump 70 and the valves 74,76, the electronic controller 88 controls the flow of liquid sterilant and rinse liquid into and out of the balloon 14 and balloon lumen 16 of each catheter 10 connected to the manifold 66 and controls the contact time of the liquid sterilant and the rinse liquid within the balloon 14 and balloon lumen 16.

The invention has been described with reference to the preferred embodiments. Obviously, modifications and alterations will occur to others upon reading and understanding the preceding detailed description. It is intended that the invention be construed as including all such modifications and alterations insofar as they come within the scope of the appended claims or the equivalents thereof.

- 22 -

Having thus described the preferred embodiments, what is claimed is:

- 1. A method of cleaning, decontaminating, and sterilizing an elongated catheter having outer surfaces, a balloon, and a dead-end balloon lumen open at a first end and closed by said balloon at a second end which is spaced from said first end, said method characterized by:
- (a) cleaning and decontaminating the outer surfaces of the catheter with a liquid decontaminant;
- (b) cleaning, decontaminating, and sterilizing the balloon lumen and the interior of the balloon with a liquid sterilant;
 - (c) drying the catheter; and,
- (d) sterilizing at least the outer surfaces of the catheter with a gaseous or plasma sterilant.
- 2. The method as set forth in claim 1 wherein the step (b) of sterilizing the balloon lumen and the interior of the balloon of the catheter is further characterized by:
- (b1) rinsing the interior of the balloon and the balloon lumen with a sterile rinse liquid to rinse any dye and saline from the balloon and the balloon lumen;
- (b2) communicating the liquid sterilant into 10 the interior of the balloon and into the balloon lumen;
 - (b3) retaining the liquid sterilant in the balloon and balloon lumen for a selected amount of time;
 - (b4) draining the liquid sterilant from the balloon and the balloon lumen; and,
- 15 (b5) repeating steps (b2) (b4) until the interior of the balloon and the balloon lumen are sterilized.

- 23 -

3. The method as set forth in claim 2 wherein the step (b5) is further characterized by repeating steps (b2) - (b4) at least eight times, and wherein step (b3) is further characterized by retaining the liquid sterilant in the balloon and in the balloon lumen for at least approximately two minutes.

5

5

time.

- 4. The method as set forth in claim 2 wherein step (b5) is further characterized by repeating steps (b2) (b4) at least ten times, and wherein step (b3) is further characterized by retaining the liquid sterilant in the balloon and in the balloon lumen for at least one and one half minutes.
- 5. The method as set forth in any of claims 1-4 wherein the catheter also includes a guide wire lumen and wherein the cleaning and decontaminating step (a) is further characterized by:
- blocking the balloon lumen to prevent the entrance of the liquid decontaminant into the balloon lumen and the interior of the balloon;

heating a bath of the liquid decontaminant to a temperature of greater than 40 degrees C;

immersing the catheter in the heated bath of liquid decontaminant for a selected amount of time; and, flushing the liquid decontaminant through the guide wire lumen of the catheter for a selected amount of

- 24 -

- 6. The method as set forth in any of claims 1-5 further characterized by repeating the cleaning and decontaminating step (a) following the sterilizing step (b).
- 7. The method as set forth in any of claims 1-6 wherein the catheter also includes a guide wire lumen and wherein the method is further characterized by, before the cleaning and decontaminating step (a):
- 5 injecting a sterile fluid into the balloon and balloon lumen of the catheter;

blocking the balloon lumen of the catheter to prevent the sterile fluid from escaping the balloon and balloon lumen; and,

- 10 treating the outer surfaces and the guide wire lumen of the catheter with an anticoagulant solution.
 - 8. The method as set forth in any of claims 1-7 wherein step (b) is further characterized by:

rinsing the balloon and balloon lumen with a sterile rinse liquid after sterilizing the balloon lumen and the interior of the balloon until the concentration of sterilant in the balloon and balloon lumen is below a predetermined threshold.

- 9. The method as set forth in any of claims 2-8 further characterized by, after step (b5):
- (b6) communicating a sterile rinse liquid into the interior of the balloon and into the balloon lumen;
- (b7) retaining the sterile rinse liquid in the balloon and balloon lumen for a selected amount of time;

5

(b8) draining the sterile rinse liquid from the balloon and the balloon lumen:

(b9) repeating steps (b6) - (b8) until the
10 amount of sterilant remaining in the balloon and the
balloon lumen is below a selected threshold.

- 10. An apparatus for sterilizing at least a dead-end balloon lumen and an interior of a balloon closing a distal end of said dead-end balloon lumen of each of a plurality of catheters using the method of any of claims 1-9, the apparatus characterized by:
 - a supply of liquid sterilant;
- a liquid bath for receiving the plurality of catheters to maintain the catheters at a selected temperature;
- a manifold in selective fluid communication with the liquid sterilant supply, the manifold including a plurality of catheter connection fittings for connecting the balloon lumen of each of the plurality of catheters in fluid communication with the manifold;
- a pump for pumping the liquid sterilant under pressure from the liquid sterilant supply into the balloon lumen and balloon of each of the plurality of catheters connected to the manifold and for evacuating liquid sterilant from the balloon and the balloon lumen of each of the plurality of catheters;

at least one valve for selectively blocking the flow of liquid sterilant to and from the manifold; and,

an electronic controller connected to the pump and the at least one valve for controlling the operation of the pump and the at least one valve, said electronic controller programmed to control the pump and the at least one valve to repeat the operations of operating the pump to fill the balloon and balloon lumen of each of the plurality of catheters connected to the manifold with

- 26 -

liquid sterilant; closing the at least one valve to retain the liquid sterilant within the balloon lumen and balloon of each of the plurality of catheters connected to the manifold for a selected amount of time; and, opening the at least one valve and operating the pump to evacuate essentially all of the liquid sterilant from within the balloon and balloon lumen of each of the plurality of catheters connected to the manifold.




Fig. 1A




Fig. 1B




Fig. 2





Fig. 4

4/7






Fig. 6

